

NEW


recomBead EBV IgG recomBead EBV IgM

Fluorescence based particle immunoassay using recombinant antigens for the detection of IgG or IgM antibodies against Epstein-Barr Virus (EBV)

The Epstein-Barr virus, an ubiquitously occurring herpes virus, can cause the symptoms of infectious mononucleosis (Pfeiffer's disease) on primary infection. Moreover, as a result of the lifelong persistence of this pathogen, reactivations can occur, especially in immuno-incompetent persons.

Due to the diversity of symptoms caused by EBV infection and their correspondence with the symptoms of other diseases, a secure EBV diagnosis is of great relevance for differential diagnostics. One of the main tasks in routine diagnosis is therefore the serological differentiation of a primary infection from a past infection and the exclusion of an EBV infection.

The modern Luminex® multiplex technology integrates the advantages from ELISA and strip assays: Quantitative detection of antibodies against individual antigens. Highly specific and characteristic EBV proteins are used for the multiplex *recomBead* EBV test systems. A classification of the different stages of infection (primary infection, past infection) is possible due to the consecutive appearance of antibodies against the different EBV antigens.




Product Advantages

- Reliable EBV screening due to recombinant antigens
 - Safe detection of past infections due to the very high specificity with the diagnostic EBV key antigen EBNA-1
 - Safe detection of acute EBV infections already in the early phase due to the optimized antigen combination in the *recomBead* EBV IgG and *recomBead* EBV IgM test systems
- Automated interpretation with feasible connection with LIMS
- Integration of advantages from ELISA and confirmation assay: Quantitative detection of antibodies against individual antigens
- Ideal screening or confirmation assay for high sample throughput
- Very high measuring accuracy and very good reproducibility of test results, therefore reliable testing of follow-up samples
- Integrated controls - no additional control samples necessary
- Small sample volume (10 µl)
- Combination of all Mikrogen *recomBead* test systems on one plate due to unified processing and exchangeable reagents
- CE label: The *recomBead* EBV test systems meet the high standard of the EC directive 98/79/EC on in vitro diagnostic medical devices

Relevant EBV antigens and recombinant analogs

EBV antigen groups	Abbreviation	Recombinant antigen	Usage
Nuclear antigen	EBNA-1	p72	IgG
Virus capsid / structural antigen	VCA	p23 p18	IgG IgG, IgM
„Immediate Early Antigen“	IEA	BZLF1	IgG, IgM
„Early Antigen“	EA	p54 p138	IgG, IgM

Test Principle and Procedure



1st Incubation

Microspheres coated with EBV specific antigens are incubated with diluted serum or plasma for **60 min.**

wash 3 times

2nd Incubation

Phycoerythrin marked anti-human antibodies (IgG or IgM specific) are added. Incubate for **30 min.**

Aspirate and add system fluid

Measurement

Either with Luminex® 100™ or Luminex® 200™ system

Evaluation

Diagnostic Sensitivity

	<i>recomBead EBV IgG</i>				<i>recomBead EBV IgM</i>
	EBNA-1	VCA		EA + IEA	
	Past infection* (n=499)	Acute infection (n=112)	Past infection (n=575)	Acute infection (n=112)	Acute infection (n=92)**
negative	0	71	1	5	0
positive	499	41	574	107	92
Sensitivity	100 %	36,6 %	99,8 %	95,5 %	100 %

* Patients with past EBV infection and anti-EBNA-1 IgG antibodies

** Patients with acute EBV infection and anti-EBV IgM antibodies

Diagnostic Specificity

	<i>recomBead EBV IgG</i>				<i>recomBead EBV IgM</i>
	EBNA-1		VCA	EA + IEA	
	EBV sero negative (n=142)	Acute infection (n=112)	EBV sero negative (n=142)	EBV sero negative (n=142)	EBV sero negative (n=142)
negative	142	112	142	142	142
positive	0	0	0	0	0
Specificity	100 %	100 %	100 %	100 %	100 %

Article-No

4552 ***recomBead EBV IgG***
Reagents for 96 determinations

4553 ***recomBead EBV IgM***
Reagents for 96 determinations

Storage and Shelf Life

At +2°C - +8°C
6 months from the date of production