

## Hand-Held Magnetic Plate Washer

The Hand-Held Magnetic Plate Washer is designed for use with the ProcartaPlex™ assay to simplify the bead wash process. The Hand-Held Magnetic Plate Washer provides you a low cost opinion for magnetic separation and automated washing of bead-based ProcartaPlex™ assays without loss of assay performance.

The Hand-Held Magnetic Plate Washer employs a magnet layer located around the bottom of the well capable to pull the beads to the side of the well. The magnetic frame is white polycarbonate with a corrosion resistant plate seal below, both set up to polypropylene base. The clip system provided is capable of holding a variety of microplates\*. The ring on the base plate facilitates gripping for all sizes of hands.\*



\*ProcartaPlex™ Assays are optimized to work with Greiner black Flat Bottom Plates, included in every Multiplex or Basic Kit.

\*For proper usage of Hand-Held Magnetic Plate Washer please refer to the ProcartaPlex™ Product Manual provided.

NOTRH AMERICA  
888.999.1371  
info@eBioscience.com  
tech@eBioscience.com

EUROPE  
+43 1 796 40 40-304  
europe@eBioscience.com  
tech@eBioscience.com

[www.eBioscience.com](http://www.eBioscience.com)