

# PCR Accessory Products

## Compression Pad

Package of 6 pads

An optional accessory to ensure a completely airtight system during thermal cycling.



## Adhesive PCR Seals

Package of 6 seals

On average, 5-10% extra seals are included in the kits. Customers requiring additional seals may purchase these separately.



## PCR Seal Applicator

Package of 1 applicator

Intended for easier application and sealing of Adhesive PCR Seals.



## PCR Cap Strips

Package of 120 strips

Intended for customers who prefer to use Cap Strips instead of Adhesive Seals. These caps are optimized for our flat-top PCR plates.



## PCR Plate Support Grid

Package of 6 grids

Intended for easier handling of cut PCR plates when putting these into and taking them out of the thermal cycler.



### Ordering information:

Product	Contents	Product No.
PCR Plate Support Grid	6 grids	103.501-06
Adhesive PCR Seals	6 seals	103.502-06
PCR Cap Strips	120 cap strip	103.503-120
PCR Seal Applicator	1 applicator	103.504-01
Compression Pad	6 pads	103.505-06

**OLERUP**

Distributed by Olerup  
[www.olerup.com](http://www.olerup.com)

CareDx-171205-5